



Culture flasks, Pyrex borosilicate glass

Scilabware

Culture flasks, Pyrex® borosilicate glass

- Manufactured from Pyrex® borosilicate glass
- Erlenmeyer
- Straight neck
- Suitable for metal caps



Code	Description	Packaging
LLG09141175	Culture flasks, Pyrex® borosilicate glass, 250 ml, Diam. 148 mm, Neck diam. 38 mm, Height 85 mm	10 pz
LLG09141176	Culture flasks, Pyrex® borosilicate glass, 500 ml, Diam. 183 mm, Neck diam. 38 mm, Height 105 mm	5 pz
LLG09141177	Culture flasks, Pyrex® borosilicate glass, 1000 ml, Diam. 229 mm, Neck diam. 38 mm, Height 131 mm	2 pz
LLG09141178	Culture flasks, Pyrex® borosilicate glass, 2000 ml, Diam. 302 mm, Neck diam. 38 mm, Height 166 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

