



Accesorios para la determinación de DQO

Behr

Accesorios para la determinación de DQO



Code	Description	Packaging
LLG06052520	Insert-/top-frame for CSB 12 E	1 pz.
LLG06052521	Vessel holding frame xsample rack insert	1 pz.
LLG06060786	Reaction vessel RG 1	1 pz.
LLG06060787	Reaction vessel with vessel adapters	1 pz.
LLG06060788	COD air cooler	1 pz.
LLG06086631	Magnetic agitator bar set for RG	12 pz.
LLG06089904	Insert-/top-frame for CSB 24 E	1 pz.
LLG09843765	Digital burette DIT 50for titration, 50	1 pz.
LLG09920701	Insert-/top-frame E 6/B	1 pz.
LLG09920711	Universal dispenser DOS 10	1 pz.
LLG09920712	Universal dispenser DOS 50	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

