



Cubetas para fotómetro

AQUALYTIC

Cubetas para fotómetro



Code	Description	Packaging
LLG06230635	Round cuvettes with coverØ 24mm, pack/5	5 pz.
LLG09699226	Adapter for 16mm cuvettes, grey	1 pz.
LLG09699228	Adapter for 13mm cuvettes, grey	1 pz.
LLG09699231	Set with 12 round cuvettes24mm with lid	12 pz.
LLG09699240	Set of 10 round vials 16mm	10 pz.
LLG09699243	Sealing ring for cuvetteØ 24 mm, gre/12	12 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

