



Cubetas para fotómetro

AQUALYTIC

Cubetas para fotómetro



Code	Description	Packaging
LLG06230635	Round cuvettes with coverØ 24mm, pack/5	5 pz.
LLG09699226	Adapter for 16mm cuvettes, grey	1 pz.
LLG09699228	Adapter for 13mm cuvettes, grey	1 pz.
LLG09699231	Set with 12 round cuvettes24mm with lid	12 pz.
LLG09699240	Set of 10 round vials 16mm	10 pz.
LLG09699243	Sealing ring for cuvetteØ 24 mm, gre/12	12 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

