




## Gradillas para tubos de ensayo, alambre/nylon

### Gradillas para tubos de ensayo, alambre/nylon

Alambre con revestimiento de nylon blanco. Orificios de 18 x 18 mm. Rejilla superior a 30 mm de altura. Resistente hasta 80 °C.



Code	Description	Packaging
LLG09193762	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193764	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193766	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193768	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193772	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193774	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193776	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193788	Test tube stands,wire,nylon-coated,18x1	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

