




## Gradillas para tubos de ensayo, alambre/nylon

### Gradillas para tubos de ensayo, alambre/nylon

Alambre con revestimiento de nylon blanco. Orificios de 18 x 18 mm. Rejilla superior a 30 mm de altura. Resistente hasta 80 °C.



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG09193762 | Test tube stands,wire,nylon-coated,18x1 | 1 pz.     |
| LLG09193764 | Test tube stands,wire,nylon-coated,18x1 | 1 pz.     |
| LLG09193766 | Test tube stands,wire,nylon-coated,18x1 | 1 pz.     |
| LLG09193768 | Test tube stands,wire,nylon-coated,18x1 | 1 pz.     |
| LLG09193772 | Test tube stands,wire,nylon-coated,18x1 | 1 pz.     |
| LLG09193774 | Test tube stands,wire,nylon-coated,18x1 | 1 pz.     |
| LLG09193776 | Test tube stands,wire,nylon-coated,18x1 | 1 pz.     |
| LLG09193788 | Test tube stands,wire,nylon-coated,18x1 | 1 pz.     |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5ª Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

