




Gradillas para tubos de ensayo, alambre/nylon

Gradillas para tubos de ensayo, alambre/nylon

Alambre con revestimiento de nylon blanco. Orificios de 18 x 18 mm. Rejilla superior a 30 mm de altura. Resistente hasta 80 °C.



Code	Description	Packaging
LLG09193762	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193764	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193766	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193768	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193772	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193774	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193776	Test tube stands,wire,nylon-coated,18x1	1 pz.
LLG09193788	Test tube stands,wire,nylon-coated,18x1	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5ª Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

