


**Tela metálica de alambre****RSG****Tela metálica de alambre**

Alambre de acero inoxidable liso con centro resistente al calor, malla: 0,6mm x 2,0mm



Code	Description	Packaging
LLG09033653	Retine spargifiamma Largh. 135 mm Lung	1 pz.
LLG09033655	Retine spargifiamma Largh. 155 mm Lung	1 pz.
LLG09033657	Retine spargifiamma Largh. 175 mm Lung	1 pz.
LLG09033659	Retine spargifiamma Largh. 195 mm Lung	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6ª Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

