




**Tela metálica de alambre**

**RSG**

**Tela metálica de alambre**

Alambre de acero inoxidable liso con centro resistente al calor, malla: 0,6mm x 2,0mm



Code	Description	Packaging
LLG09033653	Retine spargifiamma Largh. 135 mm Lung	1 pz.
LLG09033655	Retine spargifiamma Largh. 155 mm Lung	1 pz.
LLG09033657	Retine spargifiamma Largh. 175 mm Lung	1 pz.
LLG09033659	Retine spargifiamma Largh. 195 mm Lung	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

