



Test tube rack

Heathrow Scientific



Code	Description	Packaging
LLG06243203	Centrifuge tube rackfor 15x15ml tubes12	1 pz.
LLG09194021	Test tube rack 2x6wire	1 pz.
LLG09194022	Test tube rack 5x10wire	1 pz.
LLG09194023	Test tube rack 2x6wire	1 pz.
LLG09194024	Test tube rack 5x10wire	1 pz.
LLG09194025	Centrifuge bottle rack 250 ml	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

