



Vasos de molienda para CryoMill

Retsch

Vasos de molienda para CryoMill



Code	Description	Packaging
LLG09739287	Grinding beaker 5 mlhardened steel,for	1 pz.
LLG09739609	Grinding jar 10 mlstainless steel	1 pz.
LLG09739610	Grinding jar 25 mlzirconia	1 pz.
LLG09739611	Grinding jar 25 mlstainless steel	1 pz.
LLG09739612	Grinding jar 35 mlstainless steel	1 pz.
LLG09739613	Grinding jar 50 mlstainless steel	1 pz.
LLG09739614	Grinding jar 25 mlhardened steel	1 pz.
LLG09739615	Grinding jar 35 mlhardened steel	1 pz.
LLG09739616	Grinding jar 50 mlhardened steel	1 pz.
LLG09739617	Grinding jar 25 mlPTFE	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

