



Vasos de precipitados, acero inoxidable

Bochem

Vasos de precipitados, acero inoxidable

Especialidad: Fondo de contacto para placas calefactoras y agitadores magnéticos.



Code	Description	Packaging
LLG09013610	Beaker 100 ml, 18/10 steelwith rim and	1 pz.
LLG09013611	Beaker 250 ml, 18/10 steelwith rim and	1 pz.
LLG09013612	Beaker 500 ml, 18/10 steelwith rim and	1 pz.
LLG09013613	Beaker 1000 ml, 18/10 steelwith rim and	1 pz.
LLG09013614	Beaker 2000 ml, 18/10 steelwith rim and	1 pz.
LLG09013615	Beaker 3000 ml, 18/10 steelwith rim and	1 pz.
LLG09013616	Beaker 5000 ml, 18/10 steelwith rim and	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

