



**Cucharas dosificadoras LLG, HDPE**

**LLG Labware**

**Cucharas dosificadoras LLG, HDPE**

Para requisitos industriales elevados. Adecuadas para alimentos.



Code	Description	Packaging
LLG06240471	LLG-Measuring scoop1 liter, HDPE	1 pz
LLG09201091	LLG-Measuring scoop250 ml, HDPE	1 pz
LLG09201092	LLG-Measuring scoop500 ml, HDPE	1 pz
LLG09201093	LLG-Measuring scoop1500 ml, HDPE	1 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

