



Crystallising dishes, flat bottom, Pyrex

Crystallising dishes, flat bottom, Pyrex®

- Complies to DIN 12337
- Manufactured from Pyrex® borosilicate glass
- Suitable for use up to 400°C
- Chemically and thermally resistant



Code	Description	Packaging
LLG09000074	Crystallising dishes, flat bottom, Pyrex®, Without spout, 100 ml, Diam. 70 mm, Height 40 mm	10 pz
LLG09000075	Crystallising dishes, flat bottom, Pyrex®, Without spout, 500 ml, Diam. 115 mm, Height 65 mm	10 pz
LLG09000076	Crystallising dishes, flat bottom, Pyrex®, Without spout, 900 ml, Diam. 140 mm, Height 75 mm	10 pz
LLG09000077	Crystallising dishes, flat bottom, Pyrex®, Without spout, 2000 ml, Diam. 190 mm, Height 90 mm	10 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

