



## Estándar de referencia Pt-Co

### Estándar de referencia Pt-Co

Botellas de 500ml, certificadas de acuerdo al sistema QS, ISO 9001. Vida útil: 12 meses. Otras versiones disponibles a petición.



Code	Description	Packaging
LLG06207177	Platinum-Cobalt Colour 50 Certified C/1	1 pz.
LLG06223354	Platinum-Cobalt Colour 30 Certified C/1	1 pz.
LLG06232070	Platinum-Cobalt Colour 5 colour refer/1	1 pz.
LLG06232071	Platinum-Cobalt Colour 15 colour refe/1	1 pz.
LLG06240666	Platinum-Cobalt Colour 10 colour refe/1	1 pz.
LLG06240667	Platinum-Cobalt Colour 100 colour ref/1	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

