

**Conversion Kit for HACH* Photometer****Macherey Nagel**

With the 1" adapter for the DR3900 and DR6000 spectrophotometers, NANOCOLOR tube tests can be easily evaluated on HACH* photometers. Using the supplied programming data, the NANOCOLOR methods can be conveniently imported and the evaluation of the NANOCOLOR tube tests can be started immediately.

Remark *HACH is a registered trademark of HACH company, Loveland, Colorado. The product is neither endorsed, nor sponsored by or affiliated with the HACH Company.



Code	Description	Packaging
H2O919110	Conversion Kit for HACH* Photometer	



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

