

**Conversion Kit for HACH\* Photometer****Macherey Nagel**

With the 1" adapter for the DR3900 and DR6000 spectrophotometers, NANOCOLOR tube tests can be easily evaluated on HACH\* photometers. Using the supplied programming data, the NANOCOLOR methods can be conveniently imported and the evaluation of the NANOCOLOR tube tests can be started immediately.

Remark \*HACH is a registered trademark of HACH company, Loveland, Colorado. The product is neither endorsed, nor sponsored by or affiliated with the HACH Company.



Code	Description	Packaging
H2O919110	Conversion Kit for HACH* Photometer	



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

