



Dichloromethane

Synonyms:

CH₂Cl₂**Classification****Danger**

Molecular weight

84.93 ONU**1593**

CAS

75-09-2 Transport hazard
class:**6.1**

EEC-N:

200-838-9 Packing**III****Dichloromethane****> RS - ERBAdry - Anhydrous solvents -
Stabilized with ethanol****Special Grade**

Refractive index at 20°C	1.422 - 1.426	-
Water (K.F.)	<= 50	ppm
Residue on evaporation	<= 10	ppm
Colour	<= 10	APHA
Assay (GC) (without stabilizer)	>= 99.95	%
Ethanol	0.1 - 0.4	%
Acidity (HCl)	<= 3	ppm



Dichloromethane		> RS - ERBAdry - Anhydrous solvents - Stabilized with ethanol	Special Grade
Carbon tetrachloride	<= 20	ppm	
Chloroform	<= 50	ppm	

ID	Size	Packaging
P02910ED03	100 ml	Glass bottle with septum
P02910ED16	1 l	Glass bottle with septum