



n-Hexane 99%

Synonyms:

CH₃(CH₂)₄CH₃

Classification

Danger

Molecular weight

86 ONU

1208



CAS

110-54-3 Transport hazard class:

3

EEC-N:

203-777-6 Packing

II

n-Hexane 99%

> RS - SPECTROSOL - For optical spectroscopy

Special Grade

Clear, colourless liq. appearance	Conform	-
Identification	Conform	-
Colour	<= 10	APHA
Refractive index at 20°C	1.373 - 1.377	-
Water (K.F.)	<= 150	ppm
Total sulphur (S)	<= 5	ppm
Aromatic compounds	<= 5	ppm



Residue on evaporation	<= 5	ppm
Assay (GC)	>= 99	%
UV Transmittance	-	-
at 210 nm	>= 60	%
at 215 nm	>= 70	%
at 220 nm	>= 80	%
at 230 nm	>= 94	%
at 245 nm	>= 98	%

ID	Size	Packaging
000000000000447051	1 l	Glass bottle
000000000000447052	2.5 l	Glass bottle