



## Multiparameter meters MultiLine 3510 IDS

### Multiparameter meters MultiLine® 3510 IDS

Professional digital multiparameter meters for mobile measurement, with 3, 2 or 1 universal measurement channels, wireless preparation, colour graphic display (black/white display for Multi 3510), data logger manually or automatically (timed) and 2 (1 for Multi 3510) waterproof USB interfaces suitable for field use. Easy to clean with a soft brush under water jet.

**Multi 3630:** The Multi 3630 IDS with three universal measurement channels is suitable for all WTW IDS pH, ORP, dissolved oxygen sensors, conductivity cells and turbidity sensors, whether connected to cable or wireless. Documentation of meter and sensor data. User control over MultiLab® users. With USB-A and Mini USB-B interface and colour display.

**Multi 3620:** The Multi 3620 IDS with two universal measurement channels is suitable for all WTW IDS pH, ORP, dissolved oxygen sensors, conductivity cells and turbidity sensors, whether connected to a cable or wireless. Documentation of meter and sensor data. User control over MultiLab® users. With USB-A and Mini USB-B interface and colour display.

**Multi 3510:** Sequential multi-parameter measurement. All WTW IDS pH, ORP, dissolved oxygen sensors, conductivity cells and turbidity sensors can be connected to the universal channel, either with cable or wireless module. Data transfer to PC with MultiLab® Importer. With Mini USB-B interface and black/white display.

**Multi 3510 IDS:** Single meter in a carrying case with buffer 4, 7, KCl 3 mol / l, conductivity standard, batteries, USB cable

**Multi 3510 SET 1:** Meter in a carrying case with pH electrode SenTix® 940, buffer 4, 7, KCl 3 mol / l, conductivity standard, batteries, USB cable

**Multi 3510 SET 3:** Meter in a carrying case with conductivity sensor TetraCon® 925, conductivity standard, batteries, USB cable

**Multi 3510 IDS SET 4:** Meter in a carrying case with optical oxygen sensor FDO® 925, batteries, USB cable

**Multi 3510 SET KS1:** Single meter in a carrying case with buffer 4, 7, KCl 3 mol / l, conductivity standard, batteries, USB cable, arming SM Pro



Code	Description	Packaging
LLG06282979	Multiparameter meters MultiLine® 3510 IDS, Type Multi 3510, 1 channel, meter only	1 pz.
LLG06282978	Multiparameter meters MultiLine® 3510 IDS, Type Multi 3510 SET 1, 1 channel, SenTix® 940	1 pz
LLG06283710	Multiparameter meters MultiLine® 3510 IDS, Type Multi 3510 SET 3, 1 channel, TetraCon® 925	1 pz
LLG06270393	Multiparameter meters MultiLine® 3510 IDS, Type Multi 3510 SET 4, 1 channel, with FDO® 925	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

