

**Digital vertical autoclaves, AES series****Raypa****Digital vertical autoclaves, AES series**

These autoclaves are microprocessor controlled and monitored with 10 sterilisation programmes of which 4 are fixed and 6 are freely programmable.

- Operating temperature range 100 °C to 134 °C
- On request also available with vacuum pump and steam generator
- Agar modus (40 °C to 80 °C)
- RS 232 interface
- Timer
- Programmable auto-start (24 h)
- Fully automatic functions
- Case and lid shrouding are AISI 304 stainless steel
- Sterilisation chamber, lid and sealing system are AISI 316 stainless steel
- Electrical heating with protected heating elements in the interior of the chamber
- Safety pressure valve
- Safety thermostat
- Automatic steam controlled by microprocessor
- Manual valves for emptying the sterilisation chamber and for draining steam

-Supply requirements:

AES-12 and AES-75: 230 V 16 A single phase

AES-110 and AES-150: 400 V + Neutral (3-phase)



Code	Description	Packaging
LLG06267578	Digital vertical autoclaves, AES-50, 50 Litres	1 pz
LLG09851893	Digital vertical autoclaves, AES-12, 12 Litres	1 pz
LLG09851894	Digital vertical autoclaves, AES-28, 28 Litres	1 pz
LLG09851895	Digital vertical autoclaves, AES-75, 75 Litres	1 pz
LLG09851896	Digital vertical autoclaves, AES-110, 110 Litres	1 pz
LLG09851897	Digital vertical autoclaves, AES-150, 150 Litres	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

