



## Sample bottles, PP, for water sampling, sterile

### Sample bottles, PP, for water sampling, sterile

These gamma sterilised bottles are specially designed for taking, handling, transporting and storing water samples. Manufactured either as clear or amber colour from polypropylene. Clear bottles are manufactured from glass-like transparent special grade polypropylene. Amber bottles are manufactured from an amber coloured clear polypropylene which is highly translucent but still preventing the UV light penetrating inside the bottle. Square body form offers easy handling, and optimum space saving solutions during storage period. Wide bottle neck (GL 50) provides easy filling and pouring. Cap is made from PE with self sealing lock with integral security ring. Unique cap inner design ensures hermetical closure and guarantees an excellent leak and contamination proof bottle. Embossed graduation ( precision of  $\pm 2\%$  ) on side wall enables volume verifications. Each bottle is labelled specifying e.g. sample and source descriptions, code and lot numbers and expiry date.

Bottles with liquid sodium thiosulfate (20 mg/l - acc to ISO/FIS 19458:2006 ) are ideal to analyse water meant for human consumption, swimming pools, and any other water where chlorine may modify the composition of the sample.

Empty bottles, without sodium thiosulfate ; ideal to sample water for physical and chemical analysis, as well as for other sample requirements that demands sterile bottles.

Single packed bottles available on request.



| Code        | Description  | Packaging |
|-------------|--|-----------|
| LLG09001077 | Sample bottles, PP, for water sampling, sterile, clear, sterile, without sodium thiosulfate, 1000 ml | 30 pz.    |
| LLG09001083 | Sample bottles, PP, for water sampling, sterile, clear, sterile, with sodium thiosulfate, 1000 ml    | 30 pz.    |
| LLG09001075 | Sample bottles, PP, for water sampling, sterile, clear, sterile, without sodium thiosulfate, 250 ml  | 108 pz.   |
| LLG09001081 | Sample bottles, PP, for water sampling, sterile, clear, sterile, with sodium thiosulfate, 250 ml     | 108 pz.   |
| LLG09001076 | Sample bottles, PP, for water sampling, sterile, clear, sterile, without sodium thiosulfate, 500 ml  | 84 pz.    |
| LLG09001082 | Sample bottles, PP, for water sampling, sterile, clear, sterile, with sodium thiosulfate, 500 ml     | 84 pz.    |
| LLG09001085 | Sample bottles, PP, for water sampling, sterile, amber, sterile, with sodium thiosulfate, 500 ml     | 84 pz.    |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

