

**Incubators INplus and IFplus****Memmert****Incubators INplus and IFplus**

Incubators for research, medicine, pharmaceuticals and food technology. Organic chamber loads require gentle heating. For this reason, the heating and control system are especially optimised for low temperatures of up to +80 °C. To prevent temperature overshoots, temperature is increased within a very narrow control range and kept exactly at the setpoint value.

- with all-round surface heating
- fresh air is pre-heated
- with inner glass door

TwinDISPLAY: ControlCOCKPIT with two colour displays

You can find a detailed description of the "TwinDISPLAY" model at www.memmert.com



Code	Description	Packaging
LLG09537971	Incubator IN75plus +20...+80°C, 74 ltr., natural air circulation, with TwinDISPLAY	1 pz.
LLG09537972	Incubator IN110plus +20...+80°C, 108 ltr., natural air circulation, with TwinDISPLAY	1 pz.
LLG09537988	Incubator IF110plus +20...+80°C, 108 ltr., laboured air circulation, with TwinDISPLAY	1 pz.
LLG09537976	Incubator IN750plus +20...+80°C, 749 ltr., natural air circulation, with TwinDISPLAY	1 pz.
LLG09537985	Incubator IF30plus +20...+80°C, 32 ltr., laboured air circulation, with TwinDISPLAY	1 pz.
LLG09537970	Incubator IN55plus +20...+80°C, 53 ltr., natural air circulation, with TwinDISPLAY	1 pz.
LLG09537973	Incubator IN160plus +20...+80°C, 161 ltr., natural air circulation, with TwinDISPLAY	1 pz.
LLG09537989	Incubator IF160plus +20...+80°C, 161 ltr., laboured air circulation, with TwinDISPLAY	1 pz.
LLG09537974	Incubator IN260plus +20...+80°C, 256 ltr., natural air circulation, with TwinDISPLAY	1 pz.
LLG09537975	Incubator IN450plus +20...+80°C, 449 ltr., natural air circulation, with TwinDISPLAY	1 pz.
LLG09537991	Incubator IF450plus +20...+80°C, 449 ltr., laboured air circulation, with TwinDISPLAY	1 pz.
LLG09537992	Incubator IF750plus +20...+80°C, 749 ltr., laboured air circulation, with TwinDISPLAY	1 pz.
LLG09537987	Incubator IF75plus +20...+80°C, 74 ltr., laboured air circulation, with TwinDISPLAY	1 pz.
LLG09537990	Incubator IF260plus +20...+80°C, 256 ltr., laboured air circulation, with TwinDISPLAY	1 pz.
LLG09537986	Incubator IF55plus +20...+80°C, 53 ltr., laboured air circulation, with TwinDISPLAY	1 pz.
LLG09537969	Incubator IN30plus +20...+80°C, 32 ltr., natural air circulation, with TwinDISPLAY	1 pz.



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

