

**Universal Ovens UNplus and UFplus****Memmert****Universal Ovens UNplus and UFplus**

The all-round genius among heating ovens covers a multitude of applications, ideally at temperatures above +50 °C. The fresh air is preheated in the pre-heating chamber to avoid temperature fluctuations.

TwinDISPLAY: ControlCOCKPIT with two colour displays

You can find a detailed description of the "TwinDISPLAY" model at www.memmert.com



Code	Description	Packaging
LLG09869694	Universal cabinet UN260plus +20...+300°C, 256 ltr. natural air circulation, with TwinDISPLAY	1 pz.
LLG09869711	Universal cabinet UF450plus +20...+300°C, 449 ltr. laboured air circulation, with TwinDISPLAY	1 pz.
LLG09869705	Universal cabinet UF30plus +20...+300°C, 32 Ltr. laboured air circulation, with TwinDISPLAY	1 pz.
LLG09869713	Universal cabinet UF1060plus +20...+300°C, 1060 ltr. forced circulation air, with TwinDISPLAY	1 pz.
LLG09869707	Universal cabinet UF75plus +20...+300°C, 74 ltr. laboured air circulation, with TwinDISPLAY	1 pz.
LLG09869692	Universal cabinet UN110plus +20...+300°C, 108 ltr. natural air circulation, with TwinDISPLAY	1 pz.
LLG09869693	Universal cabinet UN160plus +20...+300°C, 161 ltr. natural air circulation, with TwinDISPLAY	1 pz.
LLG09869708	Universal cabinet UF110plus +20...+300°C, 108 ltr. laboured air circulation, with TwinDISPLAY	1 pz.
LLG09869709	Universal cabinet UF160plus +20...+300°C, 161 ltr. laboured air circulation, with TwinDISPLAY	1 pz.
LLG09869710	Universal cabinet UF260plus +20...+300°C, 256 ltr. laboured air circulation, with TwinDISPLAY	1 pz.
LLG09869695	Universal cabinet UN450plus +20...+300°C, 449 ltr. natural air circulation, with TwinDISPLAY	1 pz.
LLG09869712	Universal cabinet UF750plus +20...+300°C, 749 ltr. laboured air circulation, with TwinDISPLAY	1 pz.
LLG09869706	Universal cabinet UF55plus +20...+300°C, 53 ltr. laboured air circulation, with TwinDISPLAY	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
 Reagents operates with
 a Certified Quality
 Management System



