

**Stirrer bearings with cooling jacket, DURAN tubing****Lenz****Stirrer bearings with cooling jacket, DURAN® tubing**

Made of DURAN® tubing. Stirrer bearings with NS cone and cooling water connections (thread GL14), bearing surface ground and polished.



Code	Description	Packaging
LLG06234956	Stirrer bearings with cooling jacket, DURAN® tubing, For diam. stirring rod 10 mm, Cone 29/32 NS	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

