



Columns acc. to Hempel with fixed vacuum jacket, DURAN tubing

Lenz

Columns acc. to Hempel with fixed vacuum jacket, DURAN® tubing

Made of DURAN® tubing. For packings like Raschig rings or glass beads, with NS socket and NS cone, with fixed vacuum jacket.



Code	Description	Packaging
LLG06241715	Columns acc. to Hempel with fixed vacuum jacket, DURAN® tubing, Fill height 200 mm, Cone 14/23 NS, Socket 14/23 NS	1 pz.
LLG09012327	Columns acc. to Hempel with fixed vacuum jacket, DURAN® tubing, Fill height 300 mm, Cone 29/32 NS, Socket 29/32 NS	1 pz.
LLG09012328	Columns acc. to Hempel with fixed vacuum jacket, DURAN® tubing, Fill height 500 mm, Cone 29/32 NS, Socket 29/32 NS	1 pz.
LLG09012329	Columns acc. to Hempel with fixed vacuum jacket, DURAN® tubing, Fill height 800 mm, Cone 29/32 NS, Socket 29/32 NS	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

