



Accessories for Universal Ovens and Incubators

Memmert

Accessories for Universal Ovens and Incubators

Additional stainless steel grids and drip trays are available on request.



Code	Description	Packaging
LLG06266337	Accessories for Universal Ovens and Incubators, Type Stainless steel grid, For Volume 1060 l, Load max. 60 kg	1 pz.
LLG06232433	Accessories for Universal Ovens and Incubators, Type Stainless steel grid, For Volume 110/160 l, Load max. 20 kg	1 pz.
LLG09537204	Accessories for Universal Ovens and Incubators, Type Stainless steel grid, For Volume 260 l, Load max. 20 kg	1 pz.
LLG09537203	Accessories for Universal Ovens and Incubators, Type Stainless steel grid, For Volume 30 l, Load max. 20 kg	1 pz.
LLG06222688	Accessories for Universal Ovens and Incubators, Type Stainless steel grid, For Volume 450/750 l, Load max. 30 kg	1 pz.
LLG06231512	Accessories for Universal Ovens and Incubators, Type Stainless steel grid, For Volume 55/75 l, Load max. 20 kg	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

