

**LLG-Magnetic Stirring Bar Retrievers, PTFE****LLG Labware****LLG-Magnetic Stirring Bar Retrievers, PTFE**

PTFE-encapsulated stirring bar retriever with strong permanent magnet (Alnico 5), universal chemical resistance. For the removal of stirring bars from aggressive liquids, prevents loss of stirring bars. Lower end diameter 10mm.



Code	Description	Packaging
LLG04008501	Magnetic Stirring Bar Retrievers, PTFE, Diam. 8 mm, Length 300 mm	1 pz
LLG09197717	Magnetic Stirring Bar Retrievers, PTFE, Diam. 8 mm, Length 150 mm	1 pz
LLG09197725	Magnetic Stirring Bar Retrievers, PTFE, Diam. 8 mm, Length 250 mm	1 pz
LLG09197735	Magnetic Stirring Bar Retrievers, PTFE, Diam. 8 mm, Length 350 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

