

**Light Microscopes Lab-Line OBF****Kern & Sohn****Light Microscopes Lab-Line OBF**

The high-performance compound microscope for every laboratory with fixed, pre-centred Koehler illumination.

- Binocular Tube (OBF 121, OBF 122, OBF 123)
- Trinocular Tube (OBF 131, OBF 132, OBF 133)
- Tube 30° inclined , 360° rotatable
- Widefield eyepiece: WF 10x/18mm
- Diopter adjustment at one eyepiece tube
- Quadruple nosepiece (backward inclined)
- Objective 4x, 10x, 40x, 100x (Oil immersion) (different Quality objectives in different series)
- Achromatic-system and achromatic objective
- Coaxial pre-and fine focussing
- ABBE 1,25 N.A. Condensor with aperture diaphragm , focusable and pre-centered koehler illumination
- Lightfield diaphragm
- Available with 20W Halogen or 3W LED illumination
- Individuell Phase Contrast Units
- Dust Cover, eyepiece rubber caps and manual instruction included
- Additional eyepiece and filters available as accessoires





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

