

**Multi-position magnetic stirrers MIXdrive with external control use in waterbath and oven**

Without heating. With non-mechanical, electronic inductive drive and dust / water tight housing, 100 % wear-free and maintenance-free. Ideal for batch series stirring and for use in water baths. Choice of 6-place stirring position models which can stir up to 6 x 1 L erlenmeyer flasks, 15-place models for up to 15 x 250 ml beakers (high form) and 60-place models for small sample containers and using appropriate test tube stands, up to 60 test tubes or centrifuge tubes. Control units available as price-optimized eco 10 W, high-power 20 W or ultra high-power 40W version. Preselectable mixing power. SoftStart. 3 years warranty. Made in Germany.

**Please order control unit separately.**

**Specifications**

**Operating conditions:** -10 °C to +50 °C at 100 % R.H. max.

**Operating temp. range**

**Standard:**

**HT-Version:**

in air and water +50 °C

in air (ovens) +200 °C, in water +95 °C

**Stirrer unit housing:** Stainless steel

**IP code:** IP 68



Code	Description	Packaging
LLG09645961	Magnetic stirring system MIXdrive 1 eco external control -Please order control unit separately-	1 pz.
LLG09645967	Magnetic stirring system MIXdrive 60 external control -Please order control unit separately-	1 pz.
LLG09645968	Magnetic stirring system MIXdrive 60 HT (+200°C), external control - Please order control unit separately-	1 pz.
LLG09645985	Magnetic stirr system MIXdrive 1 XS external controll	1 pz.
LLG09645960	Magnetic stirring system MIXdrive 1 HT (+200°C), external control - Please order control unit separately-	1 pz.
LLG09645962	Magnetic stirring system MIXdrive 1 eco HT (+200°C), external control -Please order control unit separately-	1 pz.
LLG09645963	Magnetic stirring system MIXdrive 6 external control -Please order control unit separately-	1 pz.
LLG09645965	Magnetic stirring system MIXdrive 15 external control -Please order control unit separately-	1 pz.
LLG09645986	Magnetic stirr system MIXdrive 1 XS HT (+200°C), external controll	1 pz.
LLG09645959	Magnetic stirring system MIXdrive 1 external control -Please order control unit separately-	1 pz.
LLG09645966	Magnetic stirring system MIXdrive 15 HT (+200°C), external control - Please order control unit separately-	1 pz.
LLG09645964	Magnetic stirring system MIXdrive 6 HT (+200°C), external control - Please order control unit separately-	1 pz.





CARLO ERBA Reagents S.A.S.  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



*CARLO ERBA  
Reagents operates with  
a Certified Quality  
Management System*

