

**Multi-position magnetic stirrers MIX 6 / 15 / 8 XL / 12 XL****43953****Multi-position magnetic stirrers MIX 6 / 15 / 8 XL / 12 XL**

Without heating. With non-mechanical, electronic induction drive, 100 % wear-free and maintenance-free. Ideal for batch series stirring. Choice of 6-place, stirring position models which can stir up to 6 x 1000 ml. Erlenmeyer flasks or 15-place models for up to 15 x 250 ml (high form) or 8 / 12 x 600 ml beakers (high form). Digital stirring speed display and soft-start function. With 4 adjustable power level stirring speeds (5 W to 20 W). Easy to clean, hermetically sealed casing which can be washed under running water, ultra flat, space-saving profile and 100 % synchronous speed control in all stirring positions. SoftStart. 3 years warranty, Made in Germany.

Specifications**Max. volume stirred (H₂O): 3 l****Motor rating input: 20 W****Speed range Mix 6 / 15 :****Speed range Mix 8 XL / Mix 12 XL 100 ... 2000 min⁻¹
100 ... 1600 min⁻¹****Permissible ambient temperature: -10 ... 50 °C****Permissible humidity: 80 %****Display: digital****IP code: IP 65****Power supply: 115/230 V, 50/60 Hz**

Code	Description	Packaging
LLG09645957	Multiple Magnetic stirrer MIX 15 with internal control	1 pz.
LLG09645956	Multiple Magnetic stirrer MIX 6 with internal control	1 pz.
LLG09646026	Multiple Magnetic stirrer MIX 8 XL for 8 x 600ml beaker glasses (tall form), 100-1600 rpm	1 pz.
LLG09646027	Multiple Magnetic stirrer MIX 12 XL for 12 x 600ml beaker glasses (tall form), 100-1600 rpm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

