


**Laboratory trays LaboPlast, PP****Bürkle****Laboratory trays LaboPlast®, PP**

White multi-purpose trays. Injection-moulded. Durable and resist flexing, with strong, rounded rim. Stackable and hygienic. Surface is easy-to-clean. Resistant to many acids and alkaline solutions. Food-safe. Temperature resistant up to 60 °C.



Code	Description	Packaging
LLG09200204	Laboratory trays LaboPlast®, PP, Internal dim. 180 x 240 mm, Width 310 mm, Depth 250 mm, Height 65 mm	1 pz.
LLG09200207	Laboratory trays LaboPlast®, PP, Internal dim. 400 x 500 mm, Width 634 mm, Depth 534 mm, Height 140 mm	1 pz.
LLG09200203	Laboratory trays LaboPlast®, PP, Internal dim. 130 x 180 mm, Width 230 mm, Depth 180 mm, Height 42 mm	1 pz
LLG09200205	Laboratory trays LaboPlast®, PP, Internal dim. 240 x 300 mm, Width 370 mm, Depth 310 mm, Height 75 mm	1 pz.
LLG09200206	Laboratory trays LaboPlast®, PP, Internal dim. 300 x 400 mm, Width 520 mm, Depth 420 mm, Height 120 mm	1 pz.
LLG09200208	Laboratory trays LaboPlast®, PP, Internal dim. 500 x 700 mm, Width 846 mm, Depth 648 mm, Height 160 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

