

**Waterlanding net****Windaus-Labortechnik****Waterlanding net**

A new, easy-to-handle sample taking system for manual sampling.

General purpose thanks to

- quick and easy exchange of the sample taking equipment
- simple change by pressing the spring mechanism
- extremely strong fiber glass telescopic poles, non-corrosive
- for up to 20 different sample taking devices
- four different fiber glass telescopic poles extensible from 1.5m up to 6.00m in length

Please order telescopic poles separately.

Very robust, with nylon net, for catching flees, bugs, floating particles and water insects.

Mesh width 0.8 mm. With robust aluminium ring of 200 mm diameter. Double-stitched version. With universal connecting part.



Code	Description	Packaging
LLG06263284	Waterlanding net, Type Waterlanding net, Aperture size 0.8 mm, Material nylon	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

