

**Precision balance EHA****Kern & Sohn GmbH****Precision balance EHA**

The compact all-round model for laboratory, industry and teaching.

- High precision
- Large, shock proof weighing plate made from stainless steel
- Particularly flat design
- Simple and convenient 5-key operation
- Level indicator and levelling feet as standard
- Adjusting program CAL for quick setting of the balance accuracy
- Battery operation
- Universal mains adapter EU, GB



Code	Description	Packaging
LLG04671681	Precision balance EHA 500-1, Weighing range 500 g, Readability 0.1 g	1 pz
LLG04671682	Precision balance EHA 500-2, Weighing range 500 g, Readability 0.01 g	1 pz
LLG04671683	Precision balance EHA 3000-0, Weighing range 3000 g, Readability 1 g	1 pc
LLG04671684	Precision balance EHA 3000-1 , Weighing range 3000 g, Readability 0.1 g	1 pc



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

