



Unstirred water baths SUB Aqua Pro series, digital

Grant

Unstirred water baths SUB Aqua Pro series, digital

High quality and excellent temperature stability, in a value-for-money package designed to meet the needs of the world's researchers. The SUB Aqua Pro range is composed of eight models, including shallow and dual baths.

- Three programmable temperature presets
 - Drain tap on SUB Aqua Pro 12, 18, 26 and 34
 - Front panel lock-out
 - Countdown timer with audible buzzer
 - Ambient +5 °C to 99 °C operation
 - Stability ± 0.2 °C
 - Digital PID control for quick heat-up and precision control throughout the temperature range
 - User-settable sample protection and fixed thermal cut-out
 - High grade stainless steel tanks
 - Non drip polycarbonate lid and base tray included as standard
- Available as 120V or 230V.

Accessory polypropylene spheres or a lid must be used at temperatures between +60 °C and +99 °C. **Please order PP spheres separately.**



Code	Description	Packaging
LLG09905861	Water bath SUB AQUA PRO ADVANCED SAP2S digital, 2 L, shallow amb.+5 to 99°C, includes clear lid and base tray	1 pz.
LLG09905862	Water bath SUB AQUA PRO ADVANCED SAP5 digital, 5 L, amb.+5 to 99°C, includes clear lid and base tray	1 pz.
LLG09905860	Water bath SUB AQUA PRO ADVANCED SAP2 digital, 2 L, amb.+5 to 99°C, includes clear lid and base tray	1 pz.
LLG09905863	Water bath SUB AQUA PRO ADVANCED SAP12 digital, 12 L, amb.+5 to 99°C, includes clear lid and base tray	1 pz.
LLG09905867	Water bath SUB AQUA PRO ADVANCED SAPD digital, 5 & 12 L, dual amb.+5 to 99°C, includes clear lid and base tray	1 pz.
LLG09905865	Water bath SUB AQUA PRO ADVANCED SAP26 digital, 26 L, amb.+5 to 99°C, includes clear lid and base tray	1 pz.
LLG09905864	Water bath SUB AQUA PRO ADVANCED SAP18 digital, 18 L, amb.+5 to 99°C, includes clear lid and base tray	1 pz.
LLG09905866	Water bath SUB AQUA PRO ADVANCED SAP34 digital, 34 L, amb.+5 to 99°C, includes clear lid and base tray	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

