



Digital water baths JB Nova

Grant

Digital water baths JB Nova

General purpose water baths with stable temperature control, simple controls and fast heat up. A choice of four models with a base tray and lid included as standard.

- Ambient +5 °C to 95 °C
- Unique Set and Forgetä technology, fast reliable temperature control
- Stability ± 0.5 °C
- Simple, intuitive controls, quick and easy to set temperature
- Drain tap on 12L, 18L and 26L baths
- Practical front panel lock - disables front panel controls preventing unintentional temperature changes
- 3-year warranty

Accessory polypropylene spheres or a lid must be used at temperatures between +60 °C and +99 °C. **Please order PP spheres separately.**



Code	Description	Packaging
LLG09905870	Water bath JB NOVA GENERAL PURPOSE JBN5 digital, 5 L, amb.+5 to 95°C, includes clear lid and base tray	1 pz.
LLG09905873	Water bath JB NOVA GENERAL PURPOSE JBN26 digital, 26 L, amb.+5 to 95°C, includes clear lid and base tray	1 pz.
LLG09905871	Water bath JB NOVA GENERAL PURPOSE JBN12 digital, 12 L, amb.+5 to 95°C, includes clear lid and base tray	1 pz.
LLG09905872	Water bath JB NOVA GENERAL PURPOSE JBN18 digital, 18 L, amb.+5 to 95°C, includes clear lid and base tray	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

