

**Basic digital water baths JB Academy****Grant****Basic digital water baths JB Academy**

An ideal choice for schools and colleges requiring a basic simple-to-use quality water bath. Base tray included as standard. A great value range consisting of three models.

- Ambient + 5 °C to 95 °C
- Unique Set and Forgetä technology, fast heat-up, reliable temperature control
- Stability ± 0.5 °C
- Simple, intuitive controls, quick and easy to set temperature
- Practical front panel lock, disables front panel controls preventing unintentional temperature changes
- 3 year warranty

Accessory polypropylene spheres or a lid must be used at temperatures between +60 °C and +99 °C. **Please order lid and PP spheres separately.**



Code	Description	Packaging
LLG09905875	Water bath JB ACADEMY BASIC JBA5 digital, 5 L, amb.+5 to 95°C, includes base tray	1 pz.
LLG09905876	Water bath JB ACADEMY BASIC JBA12 digital, 12 L, amb.+5 to 95°C, includes base tray	1 pz.
LLG09905877	Water bath JB ACADEMY BASIC JBA18 digital, 18 L, amb.+5 to 95°C, includes base tray	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

