



Printers for Mettler Toledo instruments

Mettler-Toledo Online

Printers for Mettler Toledo instruments**Suitable for the following Mettler Toledo instruments:**

RS-P25, RS-P26, RS-P28: SevenExcellenceă, SevenCompactă, FiveEasy Plusă, SevenMultiă

RS-P25/USB-P25: SevenExcellenceă

RS-P26, RS-P28: SevenGo proă, SevenGo pro Duoă, Densito, Refracto, SevenEasyă

Supplied with: Printer, power cable (country specific), interface cable, paper roll (1 pcs), ribbon cartridge (1 pcs)**Features For model****Built-in real-time clock, Zero/Tare key:** RS-P26, RS-P28**Statistic application, totalling application:** RS-P28**RS232 Schnittstelle:** RS-P25, RS-P26, RS-P28**Printer unit:** Dot matrix printer**Ribbon cartridge:** Exchangeable, black ink**Paper roll:** Standard paper 57.5 mm (W), max. Ø 60 mm**Print quality:** Thermally stable and lightfast paper (GxP)**Interface:** RS232C, 9pin male connector or USB* (*SevenExcellence ONLY)**Dimensions (W x D x H):** 120 x 203 x 73mm**Net weight:** 740g (paper roll included)

Code	Description	Packaging
LLG06251155	Printers for Mettler Toledo instruments, Type RS-P26	1 pz.
LLG09773710	Printers for Mettler Toledo instruments, Type RS-P28	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

