

**Insulated container Magic Touch 2****Bel-Art Products****Insulated container Magic Touch 2ä**

Superior insulating properties combined with high chemical resistance make Magic Touch 2ä laboratory ice buckets and pans compatible not only with wet ice, but also with ultra-low temperature materials including dry ice, ice-salt mixtures, dry ice solvent mixtures, or liquid nitrogen. Their lightly textured non-slip finish and raised edge provides secure transport and a stylish look.

Expanded urethane icewares are

- Lightweight
- Durable
- Non-sweating
- Impervious to moisture and odors
- Able to withstand a wide temperature range: -196°C to 100°C
- Provided with a convenient pour spout, drip-proof rims and fill line
- Furnished with perfect fitting lid
- Easy and secure to stack because of a keyed lid and base



Code	Description	Packaging
LLG09010632	Insulated container Magic Touch 2, Type Mini Ice Pan, Colour Blue, 1.0 L, dim. (l x w x h) 183 x 183 x 116 mm	1 pz.
LLG09010633	Insulated container Magic Touch 2, Type Mini Ice Pan, Colour Red, 1.0 L, dim. (l x w x h) 183 x 183 x 116 mm	1 pz.
LLG09010635	Insulated container Magic Touch 2, Type Handy Ice Bucket, Colour Red, 2.5 L, dim. (l x w x h) (diam. x H) 268 x 164 mm	1 pz.
LLG09010636	Insulated container Magic Touch 2, Type N'Icer Ice Bucket, Colour Blue, 4.0 L, dim. (l x w x h) (diam. x H) 268 x 212 mm	1 pz.
LLG09010637	Insulated container Magic Touch 2, Type N'Icer Ice Bucket, Colour Red, 4.0 L, dim. (l x w x h) (diam. x H) 268 x 212 mm	1 pz.
LLG09010638	Insulated container Magic Touch 2, Type Midi Ice Pan, Colour Blue, 4.0 L, dim. (l x w x h) 304 x 229 x 200 mm	1 pz.
LLG09010639	Insulated container Magic Touch 2, Type Midi Ice Pan, Colour Red, 4.0 L, dim. (l x w x h) 304 x 229 x 200 mm	1 pz.
LLG09010640	Insulated container Magic Touch 2, Type Maxi Ice Pan*, Colour Blue, 9.0 L, dim. (l x w x h) 420 x 303 x 200 mm	1 pz.
LLG09010634	Insulated container Magic Touch 2, Type Handy Ice Bucket, Colour Blue, 2.5 L, dim. (l x w x h) (diam. x H) 268 x 164 mm	1 pz.
LLG09010641	Insulated container Magic Touch 2, Type Maxi Ice Pan*, Colour Red, 9.0 L, dim. (l x w x h) 420 x 303 x 200 mm	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

