

**Laboratory Bottle system DURAN YOUTILITY****DWK Life Sciences****Laboratory Bottle system DURAN® YOUTILITY**

- Specially shaped gripping zones on both sides of the bottle for easier and safer handling
- With the new YOUTILITY bottle thread opening or closing the bottle is significantly faster. The thread is fully compatible with DIN GL 45 closures and other accessories
- The slimmer YOUTILITY bottle shape allows more optimal use of limited space in autoclaves and laboratory refrigerators
- A pre-defined labelling area is compatible with the dedicated DURAN® self-adhesive YOUTILITY bottle labels
- Nominal volume is shown at the top of the easy-to-read graduation scale for fast determination of the volumes
- Glass bottle body is moulded from the tried and tested DURAN® borosilicate 3.3 Type 1 neutral glass. DURAN® glass offers very good chemical resistance and high temperature resistance
- The YOUTILITY bottle has an integrated Retrace Code to allow full product traceability. The lot specific certificate of conformance is available at: [www.duran-group.com](http://www.duran-group.com)
- Bottles are available in either clear or amber coloured glass. The amber colour meets the requirements of the USP 660> and Eur. Ph. requirements for UV light protection
- Each DURAN® YOUTILITY bottle is supplied as a complete system, with a pouring ring (cyan polypropylene) and a GL 45 cap (cyan polypropylene)



Code	Description	Packaging
LLG09071675	Laboratory Bottle system DURAN® YOUTILITY, 1000 ml, Colour amber, Ext. diam. 93 mm, Height 253 mm, Neck thread 45 GL	4 pz.
LLG09071676	Laboratory Bottle system DURAN® YOUTILITY, 1000 ml, Colour clear, Ext. diam. 93 mm, Height 253 mm, Neck thread 45 GL	4 pz.
LLG09071672	Laboratory Bottle system DURAN® YOUTILITY, 125 ml, Colour amber, Ext. diam. 55 mm, Height 124 mm, Neck thread 45 GL	4 pz.
LLG09071679	Laboratory Bottle system DURAN® YOUTILITY, 125 ml, Colour clear, Ext. diam. 55 mm, Height 124 mm, Neck thread 45 GL	4 pz.
LLG09071673	Laboratory Bottle system DURAN® YOUTILITY, 250 ml, Colour amber, Ext. diam. 66 mm, Height 158 mm, Neck thread 45 GL	4 pz.
LLG09071678	Laboratory Bottle system DURAN® YOUTILITY, 250 ml, Colour clear, Ext. diam. 66 mm, Height 158 mm, Neck thread 45 GL	4 pz.
LLG09071674	Laboratory Bottle system DURAN® YOUTILITY, 500 ml, Colour amber, Ext. diam. 78 mm, Height 193 mm, Neck thread 45 GL	4 pz.
LLG09071677	Laboratory Bottle system DURAN® YOUTILITY, 500 ml, Colour clear, Ext. diam. 78 mm, Height 193 mm, Neck thread 45 GL	4 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

