



Photographic trays LaboPlast, PVC, shallow form with ribs on bottom, profile shape rounded

**Photographic trays LaboPlast®, PVC, shallow form with ribs on bottom,
profile shape rounded**

Temperature resistant from -30 up to max. 60 °C.



Code	Description	Packaging
LLG09200269	Photographic trays LaboPlast®, shallow form without ribs on bottom, profile shape rounded, Base 110 x 160, Depth 150 mm, Width 200 mm, Height 30 mm	1 pz.
LLG09200220	Photographic trays LaboPlast®, shallow form with ribs on bottom, profile shape rounded, Base 140 x 190, Depth 180 mm, Width 230 mm, Height 40 mm	1 pz.
LLG09200270	Photographic trays LaboPlast®, shallow form without ribs on bottom, profile shape rounded, Base 140 x 190, Depth 180 mm, Width 230 mm, Height 40 mm	1 pz.
LLG09200225	Photographic trays LaboPlast®, shallow form with ribs on bottom, profile shape rounded, Base 190 x 260, Depth 240 mm, Width 320 mm, Height 50 mm	1 pz.
LLG09200275	Photographic trays LaboPlast®, shallow form without ribs on bottom, profile shape rounded, Base 190 x 260, Depth 240 mm, Width 320 mm, Height 50 mm	1 pz.
LLG09200232	Photographic trays LaboPlast®, shallow form with ribs on bottom, profile shape rounded, Base 260 x 320, Depth 320 mm, Width 370 mm, Height 50 mm	1 pz.
LLG09200282	Photographic trays LaboPlast®, shallow form without ribs on bottom, profile shape rounded, Base 260 x 320, Depth 320 mm, Width 370 mm, Height 50 mm	1 pz.
LLG09200243	Photographic trays LaboPlast®, deep form with ribs on bottom, Base 310 x 410, Depth 430 mm, Width 530 mm, Height 105 mm, Form straight	1 pz.
LLG09200283	Photographic trays LaboPlast®, deep form, without ribs on bottom, Base 310 x 410, Depth 430 mm, Width 530 mm, Height 105 mm, Form straight	1 pz.
LLG09200252	Photographic trays LaboPlast®, deep form with ribs on bottom, Base 420 x 510, Depth 520 mm, Width 610 mm, Height 110 mm, Form straight	1 pz.
LLG09200292	Photographic trays LaboPlast®, deep form, without ribs on bottom, Base 420 x 510, Depth 520 mm, Width 610 mm, Height 110 mm, Form straight	1 pz.
LLG09200267	Photographic trays LaboPlast®, deep form with ribs on bottom, Base 510 x 610, Depth 630 mm, Width 730 mm, Height 110 mm, Form straight	1 pz.
LLG09200297	Photographic trays LaboPlast®, deep form, without ribs on bottom, Base 510 x 610, Depth 630 mm, Width 730 mm, Height 110 mm, Form straight	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

