

**Vegetation scraping net****Windaus-Labortechnik****Vegetation scraping net**

A new, easy-to-handle sample taking system for manual sampling.

General purpose thanks to

- quick and easy exchange of the sample taking equipment
- simple change by pressing the spring mechanism
- extremely strong fiber glass telescopic poles, non-corrosive
- for up to 20 different sample taking devices
- four different fiber glass telescopic poles extensible from 1.5m up to 6.00m in length

Please order telescopic poles separately.

With its double reinforced aluminum ring this version can be used for scraping off particles from bushes, and trees. The net is made of nylon, mesh width 0.2mm, diameter is 280mm, net depth 500mm. With universal connecting part.



Code	Description	Packaging
LLG06263287	Vegetation scraping net, Type Vegetation scraping net, Aperture size 0.2 mm, Material nylon	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

