



Dispensing scoops, stainless steel

Dispensing scoops, stainless steel

Made from stainless steel Remanit® 4301.



Code	Description	Packaging
LLG06080286	Dispensing scoops, stainless steel, 2000 ml, Diam. 150 mm, Length Scoop 250 mm, Overall length 370 mm	1 pz.
LLG09201232	Dispensing scoops, stainless steel, 1000 ml, Diam. 105 mm, Length Scoop 190 mm, Overall length 320 mm	1 pz.
LLG09201219	Dispensing scoops, stainless steel, 120 ml, Diam. 50 mm, Length Scoop 100 mm, Overall length 190 mm	1 pz.
LLG09201235	Dispensing scoops, stainless steel, 1500 ml, Diam. 120 mm, Length Scoop 240 mm, Overall length 350 mm	1 pz.
LLG09201222	Dispensing scoops, stainless steel, 450 ml, Diam. 75 mm, Length Scoop 145 mm, Overall length 250 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

