



LLG-Timer Triple pro

LLG Labware

LLG-Timer Triple pro

The **LLG-Timer Triple pro** offers **three separate timers in one instrument**. One time can be stored per channel. This Timer can be used simultaneously as real time clock and timer. The countdown timer also runs when time is visible.

- Big 6-digit LCD display
 - Digital stopwatch (counting up)
 - Short timer (counting down)
 - With memory function
 - Stopwatch and countdown up to 23 h 59 min. 59 sec.
 - Optical and acoustic alarm when timer expires
 - Incl. magnet for mounting, hanging bracket and stand (fold-out)
- Dimensions (W x H x D): 68 x 88,7 x 16,2 mm
Weight (with battery): 87 g
Power supply: 2 x AAA 1,5 V (not included)
Warranty: 12 months

Type LLG-Timer Triple pro

Measuring capacity 23:59:59 hr: min: sec

Acoustic alarm Yes

Width 66.2 mm

Depth 13.2 mm

Height 88.2 mm

Weight 0.09 kg

Resolution sec

Vibrating alert Yes

Visual alarm Yes

Accuracy \pm sec/day

Protection class

Type of auxiliary energy 2 x 1,5 V

Dimensions (W x D x H) 66,2 x 13,2 x 88,2 mm

Colour



Code	Description	Packaging
LLG06313417	LLG-Timer Triple pro (batt. not inkl.)	1 pz



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

