

**Handheld Refractometer LLG-uniREFRACTO 5 pro****LLG Labware****LLG-uniREFRACTO 5 pro**

LLG-uniREFRACTO 5 pro is a compact digital 3-scale handheld refractometer for measuring, e.g. Brix and refractive index of liquids. Two scales are pre-set to 0 - 95 % Brix and 1.33 - 1.54 RI, the third scale can be freely selected from more than 45 different scales, including Brix, refractive index, AdBlue®, °Baume, °Butyro, Ethylene Glycol %, Calcium Chloride % and specific gravity. All three scales can be freely selected and changed. The operation is quick and easy, calibration can simply be done with distilled water. The waterproof housing (IP 65) allows cleaning under running water.

- Easy switching between the three scales
- 4-digit LCD display
- Automatic temperature correction
- Sample plate made of stainless steel

Scope of supply: Refractometer, protective cover

Specifications

Resolution: 0.1 % Brix / 0.0001 RI

Accuracy: ± 0.2 % Brix / ± 0.0003 RI

Operating temperature: 5 ... 40 °C

Measuring temperature sample: 5 ... 80 °C

Rel. Humidity: < 95 %

Dimensions (W x D x H): 120 x 55 x 35 mm

Weight: 180 g

Power supply: 2 x AAA/LR03

IP code: IP 65

Warranty: 3 years



Code	Description	Packaging
LLG06263632	LLG-uniREFRACTO 5 pro Refractometer	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

