

**pH meter SevenCompact S210**

Mettler-Toledo Online

pH meter SevenCompactä S210

The SevenCompactä pH meter S210 is the perfect choice for a wide range of applications. At the same time, working with the S210 is like child's play. Starting a measurement or calibration has never been easier, requiring only a single keypress.

- Fast and accurate measurements
- User-friendly operation
- Menu guidance in 12 languages
- Clear, well arranged color display
- Many peripheral options
- GLP compliant

Supplied with:

S210-Basic: Instrument, electrode holder, protective cover, operating instructions, quick guide, declaration of conformity, test certificate

S210-Kit: As per S210-Basic, additionally with InLab® Expert Pro-ISM, guide to pH measurements, 2 buffer sachets for pH 4.01 / 7.00 / 9.21 / 10.00

S210-Bio: As per S210-Kit, but with InLab® Routine Pro-ISM cable kit instead of InLab® Expert Pro-ISM

S210-U: As per S210-Kit, but with InLab® Versatile Pro instead of InLab® Expert Pro-ISM

Specifications:

pH-range: -2 ... 20

pH-relative accuracy: ± 0.002

pH-resolution User-definable: 0.001 / 0.01 / 0.1

mV-range: -2000.0 ... 2000.0

mV-relative accuracy: ± 0.2

mV-resolution User-definable: 0.1 / 1

Temperature range

MTC: -30.0 ... 130.0 °C

ATC: -5.0 ... 130.0 °C

Temperature accuracy: ± 0.1 °C

Display: 4.3 inch, TFT color

Power supply: 9 ... 12V / 10W

IP code: IP 54



Code	Description	Packaging
LLG09704156	pH meter SevenCompact S210, Type S210-Basic	1 pz.
LLG09704158	pH meter SevenCompact S210, Type S210-Bio	1 pz.
LLG09704157	pH meter SevenCompact S210, Type S210-Kit	1 pz.
LLG09704159	pH meter SevenCompact S210, Type S210-U	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

