




Fabric-reinforced tubing, PVC

Fabric-reinforced tubing, PVC

Transparent. PVC with terylene mesh reinforcement. Certificated by German TÜV with due regard to regulatory of DIN EN ISO 5774:2008. Proofed by BGVO directive 2002/72/EC bis 2008/39/EC for food contact. Durometer 77 Shore A.



Code	Description	Packaging
LLG09205816	Fabric-reinforced tubing, PVC, Int. diam. 6 mm, Ext. diam. 12 mm, Wall thickness 3 mm	1 pz
LLG09205818	Fabric-reinforced tubing, PVC, Int. diam. 8 mm, Ext. diam. 14 mm, Wall thickness 3 mm	1 pz
LLG09205820	Fabric-reinforced tubing, PVC, Int. diam. 10 mm, Ext. diam. 16 mm, Wall thickness 3 mm	1 pz
LLG09205823	Fabric-reinforced tubing, PVC, Int. diam. 12.5 mm, Ext. diam. 18.5 mm, Wall thickness 3 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

