



## Flow indicators LaboPlast


Bürkle

**Flow indicators LaboPlast®**

Gives a visual indication of liquid or gas flow in tubing. Operates in both directions and at any angle. Can accommodate flow rates down to approx. 150 ml/min Length incl. nozzle 88 mm. For tubing from 6mm to 11 mm diameter. Maximum pressure: 2 bar. Food Safe.

**Choice of models:**

- LiquiMobil made of SAN up to 30 °C working temperature
- CheMobil made of PMP up to 60 °C working temperature, autoclavable good, chemically resistance



Code	Description	Packaging
LLG09003912	Flow indicators LaboPlast®, Type LiquiMobil/SAN, Colour black	1 pz.
LLG09003914	Flow indicators LaboPlast®, Type CheMobil/PMP, Colour black	1 pz.
LLG09003910	Flow indicators LaboPlast®, Type LiquiMobil/SAN, Colour blue	1 pz.
LLG09003913	Flow indicators LaboPlast®, Type CheMobil/ PMP, Colour blue	1 pz.
LLG09003901	Flow indicators LaboPlast®, Type CheMobil/ PMP, Colour red	1 pz.
LLG09003911	Flow indicators LaboPlast®, Type LiquiMobil/SAN, Colour red	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

