



Flow indicators LaboPlast


Bürkle

Flow indicators LaboPlast®

Gives a visual indication of liquid or gas flow in tubing. Operates in both directions and at any angle. Can accommodate flow rates down to approx. 150 ml/min Length incl. nozzle 88 mm. For tubing from 6mm to 11 mm diameter. Maximum pressure: 2 bar. Food Safe.

Choice of models:

- LiquiMobil made of SAN up to 30 °C working temperature
- CheMobil made of PMP up to 60 °C working temperature, autoclavable good, chemically resistance



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09003912 | Flow indicators LaboPlast®, Type LiquiMobil/SAN, Colour black | 1 pz |
| LLG09003914 | Flow indicators LaboPlast®, Type CheMobil/PMP, Colour black | 1 pz |
| LLG09003910 | Flow indicators LaboPlast®, Type LiquiMobil/SAN, Colour blue | 1 pz |
| LLG09003913 | Flow indicators LaboPlast®, Type CheMobil/ PMP, Colour blue | 1 pz |
| LLG09003901 | Flow indicators LaboPlast®, Type CheMobil/ PMP, Colour red | 1 pz |
| LLG09003911 | Flow indicators LaboPlast®, Type LiquiMobil/SAN, Colour red | 1 pz |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

