

**Mini centrifuge LLG-uniCFUGE 3 pro****LLG Labware****Mini centrifuge LLG-uniCFUGE 3 pro**

This compact Mini centrifuge offers versatility, simplicity, safety and ease of usage in just one package and is specially suitable for microfiltration, cell separations and spin down applications.

- Easy rotor change via magnetic fixing
- Microprocessor controlled regulation for load-independent accuracy and integrated imbalance control
- Brushless and non-wearing high-efficiency motor
- Digital display
- Integrated timer
- Starts and stops in seconds
- Electronic safety brake on lid opening for user safety
- Lid safety interlock
- Smooth and silent, discreet running

The following types of vessel can be used:

- 8 x 1.5 ml micro tubes
- 8 x 2 ml micro tubes
- 16 x 0.2 ml PCR tubes

By using the fitting adapter:

- 8 x 0.2 ml micro tubes
- 8 x 0.5 ml micro tubes



Scope of supply: Centrifuge, rotor for 8 x 1.5/2 ml tubes, rotor for 16 x 0.2 ml PCR tubes, adapter for 8 x 0.2/0.5 ml tubes, power adapter (24 V, 1,5 A)

Specifications

- Max. RCF: 2000 x g
- Max. speed: 1000 ... 6000 rpm
- Motor: brushless DC motor
- Max. load: 8 x 1.5/2 ml or 16 x 0.2 ml PCR tubes
- Timer:
 - Dimensions (W x D x H): 1 ... 25 min, continuous
- 162 x 157 x 116 mm
- Weight: 1.1 kg
- Power supply: 100 ... 240 VAC, 50/60 Hz
- Warranty: 3 years



Code	Description	Packaging
LLG06263522	Mini centrifuge LLG-uniCFUGE 3 pro EU plug	1 pz
LLG06263527	Mini centrifuge LLG-uniCFUGE 3 pro UK plug	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

